

# Technique for measuring the particle retention of liquid filters to 10 nanometers

Gary Van Schooneveld and Don Grant  
CT Associates, Inc.  
7121 Shady Oak Road, Eden Prairie, MN 55344  
gary@ctassociatesinc.com




# Introduction

- Industrial requirements have driven the particle retention capabilities of liquid microfiltration and ultrafiltration processes to near or below 10 nanometers.
- International Technology Roadmap of Semiconductor (ITRS) has defined a critical particle size of 20 nm for 2013 and projects that in 2019, the critical particle size will be 10 nm.
- Optical particle counters traditionally used for measuring retention of ultrahigh-purity water (UPW) and chemical filters are beginning to reach their lower limits and are currently not able to measure below 30-40 nanometers.
- The potential interaction between the test particles such as polystyrene latex (PSL), gold and silica used for the challenge and the filter media is not well established.
- This presentation will discuss the use of a liquid nanoparticle sizing system capable of measuring both the size and concentration of nanoparticles in ultrapure water and its application for testing high-purity liquid filters.

# Outline

- Review the testing and particle measurement apparatus
- Comparison of candidate particles used for retention testing; in particular, their particle size distributions (PSD).
- Retention of candidate test particles by a 30 nm UPW filter
- Examples of filter retention using mono- and poly-distributed silica
- Summary


# Liquid nanoparticle sizing apparatus - filtration




## Operating Principle

- Nebulizer converts the hydrosol to an aerosol.
- DMA separates particles according to size.
- CPC measures concentrations of particles of each size.

Dynamic Mobility Analyzer (DMA)



Condensation Particle Counter (CPC)




# Liquid nanoparticle sizing system capabilities and benefits


- Broad operational range:
  - Single-particle quantitation from samples with particle concentrations at  $2 \times 10^8$  -  $1.0 \times 10^{15}$  particles/mL .
  - Sizing range from 2.5 to 450 nm with 64 channels per decade resolution (scan mode).
  - Cumulative concentration detection limit  $\geq 10$  nm of  $2 \times 10^6$  particles/mL (fixed channel mode).
- Actual concentrations, not relative concentrations are measured.
- Shape of the particle size distribution is not assumed- able to resolve multi-modal distributions.
- Highly sensitive to small changes in the PSD.
- Technique is independent of the optical and density properties of the particles.
- Suitable for on-line or off-line measurement.

# Examples of PSL particle size distributions<sup>1</sup>


20 nm PSL - Multiple dilution ratios




30nm PSL Suspension Stability




50 nm PSL 3050A-3490 - Multiple dilution ratios



100nm PSL - Multiple dilutions




# Size distributions of Ludox® Colloidal Silica Particles<sup>1</sup>



® Grace Davison





# Size distributions of gold nanoparticles<sup>1</sup>



Nominal Size (nm)	Claimed size (TEM)		Measured size (LNS)	
	Mean (nm)	CV (%)	Mean (nm)	CV (%)
10	9.3	< 15	8.4	13
20	20.3	< 8	20.8	7.4
30	30.3	< 8	30.5	7.3



# NIST 30 nm gold analysis via liquid nanoparticle sizing





## NIST Report of Investigation – Reference Material 8012 – 30 nm Gold Nanoparticles

Technique	Analyte Form	Mean Particle Size (nm)		Expanded Uncertainty
Atomic Force Microscopy	dry, deposited on substrate	24.9	±	1.1
Scanning Electron microscopy	dry, deposited on substrate	26.9	±	0.1
Transmission Electron Microscopy	dry, deposited on substrate	27.6	±	2.1
Differential Mobility Analysis	dry, aerosol	28.4	±	1.1
Dynamic Light Scattering	liquid suspension			
173° scattering angle (back scattering)		28.6	±	0.9
90° scattering angle		26.5	±	3.6
Small-Angle X-Ray Scattering	liquid suspension	24.9	±	1.2


# Study #1: 30 nm filter retention testing using three types of 30 nm particles

- Three, commercial 30 nm rated 10-inch cartridge filters from one manufacturer were tested.
- The cartridges were flushed until the filtrate approached the system background concentration ( $\sim 10^6/\text{mL} > 10 \text{ nm}$ ).
- The filter was challenged with  $2\text{E}8$  particles/mL ( $\sim 6 \text{ ppb}$ ).
- Testing conducted using single channel mode (cumulative concentration  $\geq 10 \text{ nm}$ )
- The UPW face velocity throughout the test was  $0.11 \text{ cm/min}$ .


# PSDs of 30-nm particles




# Retention of PSL followed by silica (Cartridge #1)



# Retention of gold followed by silica (Cartridge #2)



# Retention of multiple particle types (Cartridge #3)






## Observations

- Capture mechanism(s) other than sieving are active with gold and PSL that are not with silica. It is desirable to limit the capture mechanism to sieving.
- Silica is a particle that is found in water systems while gold and PSL are not.
- Silica is much less expensive than gold or PSL.

## Recommendation

- Use silica for retention testing.




# Study #2: Filter capture of mono-distributed colloidal silica particles

- Two filters tested; same pore size rating and same manufacturer.
- One filter challenged with 12.6 nm mono-distributed silica particles, the other with 18.8 nm mono-distributed silica.
- Prior to the particle challenge, the filters were flushed with UPW to reduce particle shedding.
- $3 \times 10^9$ /mL particle concentration was used.
- Particle concentrations measured upstream and downstream of the filter (scan mode).
- Filter particle retention was calculated as a function of particle size.



# Filter retention as a function of particle diameter.

Filtrate and challenge PSD's at start of 3E9/mL challenge  
(Run #1)




Test method work conducted in support of developing SEMI test guidelines for measuring retention of 5-15 nm filters used with UPW.



# Filter retention as a function of particle diameter.


Filtrate and challenge PSD - HS40 with 3E9/mL challenge  
(Run #2)




Test method work conducted in support of developing SEMI test guidelines for measuring retention of 5-15 nm filters used with UPW.



Filtrate and challenge PSD's at start of 3E9/mL challenge (Run #1)




Filtrate and challenge PSD - HS40 with 3E9/mL challenge (Run #2)



# Filter retention as a function of particle diameter.

Retention of different sized silica particles - 3E9 challenge



Overall HS40 LRV - 1.3  
Overall SM30 LRV - 0.7





Test method work conducted in support of developing SEMI test guidelines for measuring retention of 5-15 nm filters used with UPW.



# Study 3: Filter capture of poly-distributed colloidal silica particles


- Filter retention test using poly-distributed silica in UPW.
- Three silica size mixture: 12, 18 and 28 nm.
- Flush filter with UPW to reduce particle shedding. Measure particle concentration during rinse.
- $5 \times 10^9 / \text{mL} \geq 10 \text{ nm}$  particle concentration was used.
- Measure particle concentrations upstream and downstream of the filter (scanning mode).
- Filter particle retention was calculated as a function of particle size and loading.

# Silica challenge particle size distribution




# Filter Test Results


## Filter rinse in UPW



## Retention Data



# Filter retention analysis



# Summary

- Instrumentation and filter test methods have been demonstrated for the measurement of particle retention in liquid filters to 10 nanometers.
- Various nanoparticles were investigated and assessed for suitability as filter retention challenge materials.
- In the case of UPW filter retention testing, silica nanoparticles are recommended.
- Retention relative to particle diameter can currently be measured down to 10 nanometers.



# Acknowledgements

- Don Grant – CT Associates
- Erik Willis and TSI



# References

- (1) Grant DC (2011). “The effect of particle composition on filter removal of sub-30nm particles from UPW,” presented at the Executive Forum at ULTRAPURE WATER Micro 2011, Portland, OR.

